

Datasheet

WDSUB1 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00151525-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human WDSUB1 protein.

Immunogen: WDSUB1 (AAH29520.1, 1 a.a. ~ 384 a.a) full-length human protein.

Sequence:

MVKLIHTLADHGDDVNCCAFSFSLLATCSLDKTIRLYSL
RDFTELPHSPLKFHTYAVHCCCFSPSGHILASCSTDGT
TVLWNTENGQMLAVMEQPSGSPVRVCQFSPDSTCLA
SGAADGTVVLWNAQSYKLYRCGSVKDGSAAACAFSP
NGSFFVTGSSCGDLTVWDDKMRCLHSEKAHDLGITC
CDFSSQPVSDGEQGLQFFRLASCGQDCQVKIWIVSFT
DILARRTEHQLKQFTEDWSEEDVSTWLCAQDLKDLVG
IFKMNNIDGKELLNLTKESLADDLKIESLGLRSKVLRLKIE
ELRTKVKSLSSGIPDEFICPITRELMKDPVIASDGYSYE
KEAMENWISKKKRTSPMTNLVLPASVLTNRTLKMAIN
RWLETHQK

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 151525

Gene Symbol: WDSUB1

Gene Alias: FLJ36175, UBOX6, WDSAM1