

Datasheet

PPP1R2P3 monoclonal antibody (M01A), clone 2G11

Catalog Number: H00153743-M01A

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant PPP1R2P3.

Clone Name: 2G11

Immunogen: PPP1R2P3 (AAH66922, 1 a.a. ~ 205 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAASTASHRPIKGILKNKTSTTSSMVASAEQPRRSVDE
ELSKKSQKWDEINILATYHPADKGYGLMKIDEPSPPYH
SMMGDDDEDACRDTETTEAMAPGILAKKLAAGLEPK
YRIQEQUESSGEEDSDLSPEEREKKRQFEMRRKLHYNE
GLNIKLARQLISKDLHDDDEDEEMLETADGESMNTEES
NQGSTPSDQQNKLRS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 153743

Gene Symbol: PPP1R2P3

Gene Alias: MGC87149