

Datasheet

IBRDC1 polyclonal antibody (A01)

Catalog Number: H00154214-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant IBRDC1.

Immunogen: IBRDC1 (NP_689766, 2 a.a. ~ 100 a.a)
partial recombinant protein with GST tag.

Sequence:

KVQLGQVEIKCPITECFEFLEETTVVYNLTHEDSIKYKY
FLELGRIDSSTKPCPQCKHFTTFKKKGHIPTPSRSESK
YKIQCPCTCQFVWCFFKCHSPWHE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 154214

Gene Symbol: RNF217

Gene Alias: C6orf172, IBRDC1, MGC26996, dJ84N20.1