

Datasheet

DKFZp761P0423 polyclonal antibody (A01)

Catalog Number: H00157285-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant DKFZp761P0423.

Immunogen: DKFZp761P0423 (XP_291277, 2 a.a. ~ 101 a.a) partial recombinant protein with GST tag.

Sequence:

HQTLCLNPESLKMSACSDFVEHIWKPGSCKNCFCLRS
DHQLVAGPPQPRAGSLPPPPRLPPRPENCRLEDEGV
NSSPYSKPTIAVKPTMMSSEASDVWTE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 157285

Gene Symbol: PRAGMIN

Gene Alias: DKFZp761P0423