

Datasheet

PRAGMIN monoclonal antibody (M02), clone 1D9

Catalog Number: H00157285-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant PRAGMIN.

Clone Name: 1D9

Immunogen: PRAGMIN (XP_291277, 2 a.a. ~ 101 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

HQTLCLNPESLKMSACSDFVEHIWKPGSCKNCFCLRS
DHQLVAGPPQPRAGSLPPPPRLPPRPENCRLEDEGV
NSSPYSKPTIAVKPTMMSSEASDVWTE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 157285

Gene Symbol: PRAGMIN

Gene Alias: DKFZp761P0423