

Datasheet

C9orf98 MaxPab mouse polyclonal antibody (B01)

RP11-143F18.1

Catalog Number: H00158067-B01**Regulation Status:** For research use only (RUO)**Product Description:** Mouse polyclonal antibody raised against a full-length human C9orf98 protein.**Immunogen:** C9orf98 (NP_689785.1, 1 a.a. ~ 479 a.a) full-length human protein.**Sequence:**

MDATIAPHRIPPEMPQYGEENHIFELMQNMLEQLLIHQ
PEDPIPFMIQHLHRDNDNVPRIVILGPPASGKTTIAMWL
CKHLNSSLTLENLILNEFSYTATEARRLYLQRKTVPSA
LLVQLIQERLAEEDCIKQGWILDGIPETREQALRIQTLGI
TPRHVIVLSAPDTVLIERNLGKRIDPQTGEIYHTTFDWP
PESEIQNRMLMVPEDISELETAQKLLEYHRNIVRVIPSYP
KILKVISADQPCVDVFYQALTYVQSNHRTNAPFTPRVL
LLGPVSGSKSLQAALLAQKYRLVNVCCGQLLKEAVAD
RTTFGELIQPFPEKEMAVPDSLLMKVLSQRDQDQCIQ
KGWVLHGVPRDLDAHLLNRLGYNPNRVFFLNVPFD
SIMERLTLRRIDPVTGERYHLMYKPPPTMEIQARLLQN
PKDAEEQVKLKMDFYRNSADLEQLYGSAILNGDQD
PYTVFEYIESGIINPLPKKIP

Host: Mouse**Reactivity:** Human**Applications:** WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Storage Buffer:** No additive**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 158067**Gene Symbol:** C9orf98**Gene Alias:** DDX31, FLJ32704, FLJ36014,