

Datasheet

TTC16 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00158248-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TTC16 protein.

Immunogen: TTC16 (AAH31281, 1 a.a. ~ 873 a.a.) full-length human protein.

Sequence:

MTDSDEDALKVDQGSPRDIPKPWWIPAPKGILQHIFGT
SHVFQSIDVKPKVTGLTVPLKVREYYSRGQQCLEQA
DWETAVLLFSRALHLDLPQLVDFYALRAEAYLQLCDFSS
AAQNLRRAYSLQQDNCKHLERLTFVLYLQGQCLFEQC
AFLDALNVFSAHAELQPEKPCFRYRCMACLLALKQH
ACLTITNELKQDATTNADVYIFRARLYNFLQKPHLCYRD
LHSALLNPKHPQARMMLQKMVAQAQARQDAGILAV
QGKQLHALQRINRAIENNPLDPSLFLFRGTMYYRLQEF
DGAVEDFLKVLDMVTEDQEDMVRQAQRQLLLTYNDF
AVHCYRQGAYQEGVLLLNKALRDEQKEGLYINRGDC
FFQLGNLAFEAADYQQALALSPQDEGANTRTGLLQEK
MGFCEQRRKQFQKAENHFSTAIRHNPQKAQYYLYRA
KSRQLLNIFGARQDVATVLLNPKQPKLSLLMTNLFP
GMSVEEVLSTQIAHLARLQLEQMVEGSLQAGSPQGIV
GMLKRHELERQKALALQHSWKQGEPLIATSEELKATP
EIPQVKPGSSEGEAEAEPEEEEEKEKEKKGEKKSELIPS
KVASLSDSYLDQTSSASSMSFRTTGTSETEMSAICQE
YRSTSATAVTFSDSLLKTQSSDSGNNREALSHGPRKI
KATQGQRQSLSKTEPTQSQRNNSKTKATIHKRNSSK
TKATQSQRNNSKTRATQGGQSSSKTEATQGQRQS
SSEIATQGPRQEPSKTKTTRSPRQRPRKVKAAARGRS
WRPSKVDATQGRSRGLLRSSKTEAFYDSNWSLSKT
EYAQQGQRSSKAEGAQGKSQGMSSSTSSKAESTWG
PSPSLSKTEVDQDLTYYESV

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product

page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 158248

Gene Symbol: TTC16

Gene Alias: FLJ32780