

Datasheet

ZNF483 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00158399-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF483 protein.

Immunogen: ZNF483 (AAH65738.1, 1 a.a. ~ 244 a.a) full-length human protein.

Sequence:

MQAVVPLNKMTAISPEPQTLASTEQNEVPRVVTSGEQ
EAILRGNAADAESFRQRFRWFCYSEVAGPRKALSQ
WELCNQWLRPDIHTKEQILELLVFEQFLTILPGEIRIWW
KSQHPESSEEVVTLIEDLTQMLEEKDPVSQDSTVSQE
ENSKEDKMVTVCNTESCESITLKDVAVNFSRGEWKK
LEPFQKELYKEVLLENLRNLEFLDFPVSKLELISQLKWV
ELPWLLLEEVSKSSRLGSI

Host: Mouse

Reactivity: Human

Applications: IF, IHC-P, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 158399

Gene Symbol: ZNF483

Gene Alias: ZKSCAN16