

Datasheet

EXDL1 MaxPab mouse polyclonal antibody (B01)

Gene Alias: MGC33637

Catalog Number: H00161829-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human EXDL1 protein.

Immunogen: EXDL1 (NP_689809, 1 a.a. ~ 514 a.a) full-length human protein.

Sequence:

MEDSEFLAYVELLDEVEQGSVRAKASSVSLHAERTW
MEKMKVEDLNVCEPASPAPAPATSLNDLKYSPEE
EEVTTYVINQFQQKFGAAILHIKKQNVLSVAAEGANVC
RHGKLCWLQVATNCRVYLFDFLLGSRAFHNGLQMLE
DKRILKVIHDCRWLSDCLSHQYGILLNNVFDQVADVL
QFSMETGGYLPNCITTLQESLIKHLQVAPKYSFLEKR
QKLIQENPEVWFIRPVSPSLKILALEATYLLPLRLALLD
EMMSDLTTLVDGYLNTYREGSADRLGGTEPTCMELPE
ELLQLKDFQKQRREKAAREYRVNAQGLLIRTVLQPKKL
VTETAGKEEKVKGFLFGKNFRIDKAPSFTSQDFHGDV
NLLKEESLNKQATNPQHLPTTEGETSEDSSNKLICTK
SKGSEDQRITQKEHFMTPKHEFQASLSLKEETEQLLM
VENKEDLKCTKQAVSMSSFPQETRVSPSDTFYPIRKT
VVSTLPPCPALEKIDSWISPFNLNP

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 161829

Gene Symbol: EXDL1