

Datasheet

LOC161931 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: ADAD2

Gene Alias: FLJ00337, TENRL

Catalog Number: H00161931-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LOC161931 protein.

Immunogen: LOC161931 (AAH33491, 1 a.a. ~ 583 a.a) full-length human protein.

Sequence:

MASASQGADDDGSRRKPRLAASLQISPQPRPWRPLP
AQAQSAWGPAPAPATYRAEGGWPQVSVLRDSGPGA
GAGVGELGAARAWENLGEQMGKAPRPVPPAGLSLP
LKDPASQAVSLLTEYAASLGIFLLFREDQPPGPCFPF
SVSAELDGVVCPAGTANSKTEAKQQAALSALCYIRSQ
LENPESPQTSSRPPLAPLSVENILTHEQRCAALVSAGF
DLLLDERSPYWACKGTVAGVILEREIPRARGHVKEIYK
LVALGTGSSCCAGWLEFSGQQLHDCHGLVIARRALLR
FLFRQLLLATQGGPKGKEQSVLAPQPGPGPPFTLKPR
VFLHLYISNTPKGAARDIYLPPTSEGGLPHSPPMRLQA
HVLGQLKPVICYVAPSLCDTHVGCLSASDKLARWAVLG
LGGALLAHLVSPLYSTSLILADSCHDPPTLSRAIHRPC
LDSVLGPCLPPPYVRTALHLFAGPPVAPSEPTDTCR
GLSLNWSLGDPGIEVVVDVATGRVKANAALGPPSRLCK
ASFLRAFHQAARAVGKPYLLALKTYEAAKAGPYQEAR
RQLSLLLDQQGLGAWPSKPLVGKFRN

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 161931