

Datasheet

RNF190 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00162333-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human RNF190 protein.

Immunogen: RNF190 (AAH35021.1, 1 a.a. ~ 657 a.a) full-length human protein.

Sequence:

MLHDARDRQKFFSDVQYLRLDMQHKVDSEYQACLRQ
EYRRDPNEKKRDQFWGQETSFERSRFSSRSSKQSS
SEEDALTEPRSSIKISAFKCDKSLPAIDQTSVKQKHST
MTVRKAEKVDPSPADQAPMVLLRKRKPNLRRFTV
SPESHSPRASGDRSRQKQWPAPKVPVPRGADQVVQ
QEGLMCNTKLKRPNQERRNLVPSSQPMTENAPDRAK
KGDPSAPSQSELHPALSQAFQEKNSPQVLSEFSGPPL
TPTTVGGPRKASFRFRDEDFYSILSLNSRRESDDTEEE
TQSEECLWVGVRSPCSPSHHKRSRFGGTSTPQAKNK
NFEENAENCRGHSSRRSEPSHGSLRISNAMEPATERP
SAGQRLSQDPGLPDRESATEKDRGGSENAKKSPLSW
DTKSEPRQEVGVNAENVWSDCISVEHRPGTHDSEGY
WKDYLNSSQNSLDYFISGRPISPRSSVNSSYNPPASF
MHSALRDDIPVDLSMSSTSVHSSDSEGNSGFHVCQPL
SPIRNRTPFASAENHNYFPVNSAHEFAVREAEDTTLTS
QPQGAPLYTDLLLNPQGNLSLVDSSSSSPSRMNSEGH
LHVSGSLQENTPFTFFAVSHFPNQNDNGSRMAASGFT
DEKETSKIKADPEKLKKLQESLLEEDSEEEGD

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 162333

Gene Symbol: MARCH10

Gene Alias: FLJ35757, MARCH-X, RNF190