

## Datasheet

### PHOSPHO1 polyclonal antibody (A01)

**Catalog Number:** H00162466-A01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a partial recombinant PHOSPHO1.

**Immunogen:** PHOSPHO1 (NP\_848595, 168 a.a. ~ 266 a.a) partial recombinant protein with GST tag.

**Sequence:**

CARCPANMCKHKVLS DYLRERAH DGVHFERLFYVGD  
GANDFCPMGLLAGGDVAFPRRGYPMHR LIQEAQKAE  
PSSFRASVVPWETAADVRLHLQQVLKS

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** 50 % glycerol

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 162466

**Gene Symbol:** PHOSPHO1

**Gene Alias:** -