

## Datasheet

### ZNF433 purified MaxPab mouse polyclonal antibody (B01P)

**Entrez GeneID:** 163059**Gene Symbol:** ZNF433**Gene Alias:** FLJ40981**Catalog Number:** H00163059-B01P**Regulation Status:** For research use only (RUO)**Product Description:** Mouse polyclonal antibody raised against a full-length human ZNF433 protein.**Immunogen:** ZNF433 (AAH47412.1, 1 a.a. ~ 638 a.a) full-length human protein.**Sequence:**

MQETFRNLASIGKKWKPNIIYVEYENLRRNLRIVGERL  
FESKEGHQHGEILTQVPDDMLKTTTGVKSCSSVYG  
EVGSAHSSLNRHIRDDTGHKAYEYQEYQKPYKCKY  
CKKPFNCLSSVQTERAHSGRKLYVCEECKGKTFISHS  
NLQRHRIMHRGDGPYKCKFCGKALMFLSLYLIHKRTH  
TGEKPYQCKQCGKAFSHSSSLRIHERHTTGEKPYKCN  
ECGKAFHSSTCLHAHKRHTTGEKPYECKQCGKAFSS  
SHSFQIHERHTTGEKPYECKECKGKAFKCPSSVRRHER  
THSRKKPYECKHCGKVLSTFQNLHGMHTGEISHK  
CKICGKAFYSPSSLQTHEKTHTGEKPYKCNQCGKAFN  
SSSSFRYHERHTTGEKPYECKQCGKAFRSASLLQTHG  
RTHTGEKPYACECKGKPFNSNFSFFQIHERMHREEKPY  
ECKGYGKTFSLPSLFHRHERHTTGGKTYECKQCGRS  
FNCSSSFRYHGRHTTGEKPYECKQCGKAFRSASQLQI  
HGRHTTGEKPYECKQCGKAFGSASHLQMHGRHTTGE  
KPYECKQCGKSFGCASRLQMHGRHTTGEKPYKCKQC  
GKAFGCPSNLRRHGRHTTGEKPYKCNQCGKVFRCS  
QLQVHGHAHCIDTP

**Host:** Mouse**Reactivity:** Human**Applications:** IF, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Storage Buffer:** In 1x PBS, pH 7.4**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.