

Datasheet

ZNF567 MaxPab mouse polyclonal antibody (B01)

Entrez GeneID: 163081**Gene Symbol:** ZNF567**Gene Alias:** MGC45586**Catalog Number:** H00163081-B01**Regulation Status:** For research use only (RUO)**Product Description:** Mouse polyclonal antibody raised against a full-length human ZNF567 protein.**Immunogen:** ZNF567 (NP_689816.2, 1 a.a. ~ 616 a.a) full-length human protein.**Sequence:**

MDVMLENYCHLISVGCHMTKPDVILKLERGEEPWTSF
AGHTCLEENWKAEDFLVKFKEHQEKYSRSVVSINHKK
LVKEKSKIYEKFTLGKNPVNSKNLPPEYDTHGRILKNV
SELIISNLNPARKRLSEYNGYGKSLSTKQETTHPEVKS
HNQSARAFSHNEVLMQYQKTETPAQSFGYNDCEKSF
LQRGGLITHSRPYKGENPSVYNKKRRATNIEKKHTCN
ECGKSFCRKSVLILHQGIHSEEKPYQCHQCGNAFRRK
SYLIDHQRTHTEKPFVCNECGKSFRLKTALTDHQR
HTGEKSYECLQCRNAFRLKSHLIRHQRTHTEKPYEC
NDCGKSFRQKTTLSLHQRIHTGEKPYICEGKSFHQ
KANLTVHQRTHTEKPYICNECGKSFSQKTTLALHEKT
HNEEKPYICSECGKSFRQKTTLVAHQRTHTGEKSYEC
PHCGKAFRMKSYLIDHHRTHTEKPYECNECGKSFSQ
KTNLNLHQRIHTGEKPYVCNECGKSFRQKATLTVHQKI
HTGQKSYECPQCGKAFSRKSyliHHQRTHTEKPYKC
SECGKCFRQKTNLIVHQRTHTEKPYVCNECGKSFSY
KRNLIVHQRTHKGENIEMQ

Host: Mouse**Reactivity:** Human**Applications:** WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Storage Buffer:** No additive**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.