

Datasheet

MGC27016 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: RBM46

Gene Alias: MGC27016

Catalog Number: H00166863-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC27016 protein.

Immunogen: MGC27016 (NP_659416, 1 a.a. ~ 533 a.a) full-length human protein.

Sequence:

MNEENIDGTNGCSKVRTGIQNEAALLALMEKTGYNMV
QENGQRKFGGPPPGWEGPPPPRGCEVFVGKIPRDM
YEDELVPVFERAGKIYEFRLMMEFSGENRGYAFVMYT
TKEEAQLAIRILNNEYIRPGKFIGVCVSLDNCRLFIGAIP
KEKKKEEILDEMKKVTEGVVDVIVYPSATDKTKNRGFA
FVEYESHRAAAMARRKLIPGTFQLWGHTIQVDWADPE
KEVDEETMQRVKVLVVRNLMISTTEETIKAEFNFKPG
AVERVKLRDYAFVHFFNREDAVAAMSV MNGKCIDGA
SIEVTLAKPVNKENTWRQHLNGQISPENSENLIVFANKE
ESHPKTLGKLPTLPARLNGQHSPSPPEVERCTYPFYP
GTKLTPISMYSLSNHFNSAVMHLDYCNKNWAPPE
YYLYSTTSQDGKVLVYKIVIPAIANGSQSYFMPDKLCT
TLEDAKELAAQFTLLHLDYFNHRSSINSLSPVSATLSSG
TPSVLPYTSTRPYSYPGYPLSPTISLANGSHVGQRLCIS
NQASFF

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 166863