

Datasheet

UGT3A2 MaxPab mouse polyclonal antibody (B01)

Gene Alias: MGC119426, MGC119429

Catalog Number: H00167127-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human UGT3A2 protein.

Immunogen: UGT3A2 (AAI03925, 1 a.a. ~ 523 a.a)
full-length human protein.

Sequence:

MAGQRVLLLVGFLLPGVLLSEAAKILTISTVGGSHYLLM
DRVSQILQDHGHNVTMLNHNKRGPFMPDFKKEEKSYQ
VISWLAPEDHQREFKKSFDFFLEETLGGRGKFENLLNV
LEYLALQCSHFLNRKDIMDSLKNENFDMVIVETFDYCP
FLIAEKLKGKPFVAILSTSFGSLEFGLPIPLSYVPVFRSLL
TDHMDFWGRVKNFLMFFSFQRRQQHMQSTFDNTIKE
HFTEGSRPVLSHLLKAEWFINSDFAFDFAFARPLLNT
VYVGGLMEKPIKPVPQDLENFIAKFGDSGFVLVTLGSM
VNTCQNPEIFKEMNNAFAHLPQGVIWKCQCQSHWPKD
VHLAANVKIVDWLPQSDLLAHPISIRLFVTHGGQNSIME
AIQHGVPVMVGIPFGDQPENMVRVEAKKFGVSIQLKKL
KAETLALKMKQIMEDKRYKSAAVAASVILRSHPLSPTQ
RLVGWIDHVLQTTGGATHLKPYPVQQPWHEQYLLDVFFV
FLLGLTLGTLWLCGKLLGMAVWWLRGARKVKET

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 167127

Gene Symbol: UGT3A2