

Datasheet

MGC39715 polyclonal antibody (A01)

Catalog Number: H00169166-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant MGC39715.

Immunogen: MGC39715 (AAH31260, 242 a.a. ~ 341 a.a) partial recombinant protein with GST tag.

Sequence:

QNLELRFQYSEDSCWQWFVIYTKQAFLLSSCLKKMIS
EKMVKLAAENTEMQIEVPEQSKSKKYHIQQSQQKDYS
SFLSRKSKIKIAKGDCVFGNIKEEDL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 169166

Gene Symbol: SNX31

Gene Alias: MGC39715