

Datasheet

C9orf96 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00169436-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human C9orf96 protein.

Immunogen: C9orf96 (AAH36504.3, 1 a.a. ~ 680 a.a) full-length human protein.

Sequence:

MLGPGSNRRRPTQGERGPGSPGEPMEKYQVLYQLN
PGALGVNLVVEEMETKVKHVIKQVECMDDHYASQALE
ELMPLLLKLRHAHISVYQELFITWNGEISSLYLCLVMEFN
ELSFQEVIEDKRKAKKIIDSEWMQNVLGQVLDALAYLH
HLDIIHRNLKPSNIILISSDHCKLQDLSSNVLMTDKAKW
NIRAEEDPFRKSWMAPEALNFSFSQKSDIWSLGCIIID
MTSCSFMDGTEAMHLRKSRLQSPGSLKAVLKTMEEK
QIPDVETFRNLLPLMLQIDPSDRITIKDVVHITFLRGSFK
SSCVSLTLHRQMVPASITDMLLEGNVASILEVMQKFSG
WPEVQLRAMKRLLKMPADQLGLPWPELVEVVVTM
ELHDRVLDVQLCACSLLLHLLGQALVHHPEAKAPCNQ
AITSTLLSALQSHPEEEPLLVMVYSLLAITTTQESESLS
ELQNAGLLEHILEHLNSSLERSDVCASGLGLLWALLLD
GIIVNKAPLEKVPDLISQVLATYPADGEMAEASCGVFW
LLSLLGCIKEQQFEQVALLQSIRLCQDRALLVNNAYR
GLASLVKVSELAFAFKVVVQEEGSGLSLIKETYQLHRD
DPEVVENVGMMLVHLASYEEILPELVSSSMKALLQEIKE
RFTSSLVSDSSAFSKPGLPPGGSPQLGCTTSGGLE

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 169436

Gene Symbol: C9orf96

Gene Alias: MGC43306, RP11-244N20.8