

Datasheet

C9orf96 monoclonal antibody (M03), clone 2B5

Catalog Number: H00169436-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant C9orf96.

Clone Name: 2B5

Immunogen: C9orf96 (NP_714921, 1 a.a. ~ 99 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MLGPGSNRRRPTQGERGPGSPGEPMEKYQVLYQLN
PGALGVNLVVEEMETKVKHVIKQVECMDDHYASQALE
ELMPLLKLRHAHISVYQELFITWNGEI

Host: Mouse

Reactivity: Human

Applications: ELISA, RNAi-Ab, S-ELISA, WB-Re,
WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 169436

Gene Symbol: C9orf96

Gene Alias: MGC43306, RP11-244N20.8