

Datasheet

METT5D1 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00196074-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human METT5D1 protein.

Immunogen: METT5D1 (AAH30997, 1 a.a. ~ 212 a.a) full-length human protein.

Sequence:

MAKLHIPVMVDEVVHCLSPQKGQIFLDMTFGSGGHTK
AILQKESDIVLYALDRDPTAYALAEHLSELYPKQIRAML
GQFSQAEALLMKAGVQPGTFDGVLMDLGCSSMQLDT
PERGFSRLRKGSLDMRMDGGRSTGTCTIYPKNIRGGEA
CQENRFSNCSGTQHLP HHQNPAACQHRCRSISSLCY
LYTERLT TAIYPYCHQDFPGSSHICEQ

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 196074

Gene Symbol: METT5D1

Gene Alias: FLJ33979