

Datasheet

LOC200383 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: LOC200383

Gene Alias: KIAA1697

Catalog Number: H00200383-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LOC200383 protein.

Immunogen: LOC200383 (AAH15442, 1 a.a. ~ 551 a.a) full-length human protein.

Sequence:

MEQTKLNTILCQTFVFCYLWSLGGNLTENYYDSFDTFI
RTQFDDNPDPARLPNSGDLWSIHMDFDTKRLDPWERII
PTFKYNRDVPFFEMLVPTTDTVRYGYLMEKLLAVKHS
VLFTGITGVGKSIAKGLLNKIQESAGYVPVYLNFSAQ
SSARTQEIESKLERKRKNILGAPGNKRIVFVDDLNP
RLDRYGSQPPIELLRQYQDFGGFYDRNKLFWKEIQDV
TIISACAPPGGGRNPVTPRFIRHFSMLCLPMPSEHSLK
QIFQAILNGFLSDFPPAVKQTASSIVEASVEIYNKMSVD
LLPTPAKSHYVFNLRDLSKCVQDKHFGIAIDLEYFLNKP
IIFGDFIKFGADKADRIYDDMPDIEKTANVLQDYLDYDYN
LTNPKEVKLVFFQDAIEHVSRIARMIRQERGNALLVGV
GGTGKQSLTRLAAHICGYKCLQIELSRGYNYSFHEDL
RKLYKMAGVEDKNMVFLFTDTQIVVEEFLEDINNILNS
GEVPNLFKDELEQVLAATRPRAKEVGISEGNRDEVG
CARVVMWLYHKNVSSPFI

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 200383