

Datasheet

ZNF584 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00201514-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF584 protein.

Immunogen: ZNF584 (NP_775819.1, 1 a.a. ~ 421 a.a) full-length human protein.

Sequence:

MAGEAEAQLDPSLQGLVMFEDVTVYFSREEWGLLNV
TQKGLYRDVMLENFALVSSLGLAPSRSPVFTQLEDDE
QSWVPSWVDVTPVSRAEARRGFGLDGLCRVEDERA
HPEHLKSYRVIQHQDTHSEGKPRRHTEHGAAFP PGSS
CGQQQEVHVAEKLKFCSDCGKVFLKAFALLDHLITHSE
ERPFRCP TGRSAFKKSAHINPRKIHTGETAHVCNECG
KAFSYP SKLRKHQKVHTGIKPFKCSDCGKTFNRKDAL
VLHQRIHTGERPYECSKCGKTF SVLSTLIRHRKVHIGE
RPYECTECGKFFKYNNSFILHQRVHTGERPFECKQCG
KGYVTRSGLYQHWKVHTGERPYECSLCGKTF TTRSY
RNRHQQFHTEERSYECTECGKAFKHSSTLLQHKKVHT
PERRQEDRAHGKVVSC

Host: Mouse

Reactivity: Human

Applications: IF, IHC-P, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 201514

Gene Symbol: ZNF584

Gene Alias: FLJ39899