

Datasheet

PTPLB purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00201562-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PTPLB protein.

Immunogen: PTPLB (NP_940684.1, 1 a.a. ~ 254 a.a) full-length human protein.

Sequence:

MAAVAATAAAKGNNGGGGGRAGAGDASGTRKKKGPG
PLATAYLVIYNVVMTAGWLVIAGLVRLAYLAKGSYHSL
YYSIEKPLKFFQTGALLEILHCAIGIVPSSVLTSTFQVMS
RVFLIWAVTHSVKEVQSEDSVLLFVIAWTITEIRYSFYT
FSLLNHLPYLIKWARYTLFIVLYPMGVSGELLTIYAALPF
VRQAGLYSISLPNKYNFSFDYYAFLILIMISYIPIFPQLYF
HMIHQRRKILSHTEEHKKFE

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 201562

Gene Symbol: PTPLB

Gene Alias: -