

Datasheet

STK32A monoclonal antibody (M05), clone 2C3

Catalog Number: H00202374-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant STK32A.

Clone Name: 2C3

Immunogen: STK32A (AAH21666, 67 a.a. ~ 166 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

NVFKELQIMQGLEHPFLVNLWYSFQDEEDMFMVVDLL
LGGDLRYHLQQNVHFKEETVKLFICELVMALDYLQNQ
RIIHRDMKPDNILLDEHDTWLSYKSH

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 202374

Gene Symbol: STK32A

Gene Alias: MGC22688, YANK1