

## Datasheet

### PLEKHK1 MaxPab mouse polyclonal antibody (B01)

**Catalog Number:** H00219790-B01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human PLEKHK1 protein.

**Immunogen:** PLEKHK1 (AAH25765, 1 a.a. ~ 163 a.a)  
full-length human protein.

**Sequence:**

MEGPSLRGPALRLAGLPTQQDCNIQEKIDLEIRMREGI  
WKLLSLSTQKDQVLHAVKNLMVCNARLMAYTSELQKL  
EEQIANQTGRCDVKFESKERTACKGKIAISDIRIPLMWK  
DSDHFSNKERSRRYAIFCLFKMGANVFDTDVVNVDKTI  
TDICFENVIL

**Host:** Mouse

**Reactivity:** Human

**Applications:** Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** No additive

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 219790

**Gene Symbol:** RTKN2

**Gene Alias:** DKFZp686J10120, PLEKHK1, bA531F24.1