

Datasheet

MGC35295 monoclonal antibody (M03), clone 1D4

Catalog Number: H00219995-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant MGC35295.

Clone Name: 1D4

Immunogen: MGC35295 (AAH31610.1, 1 a.a. ~ 147
a.a) full-length recombinant protein with GST tag. MW of
the GST tag alone is 26 KDa.

Sequence:

MVRRGHVGIFFIEGGVPFWGGACFIISGSLSVAAEKNH
TSCLVRSSLGTNLSVMAAFAGTAILLMDFGVTNRDVD
RGYLAVLTIFTVLEFFTAVIAMHFGCQAIHAQASAPVIFL
PNAFSADFNIPSPAASAPPAYDNVAYAQGVV

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 219995

Gene Symbol: MS4A15

Gene Alias: FLJ34527, MGC35295