

Datasheet

MARCH8 polyclonal antibody (A01)

Catalog Number: H00220972-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant MARCH8.

Immunogen: 40976 (NP_659458, 77 a.a. ~ 153 a.a)
partial recombinant protein with GST tag.

Sequence:

QDICRICHCEDDESPLITPCHCTGSLHFVHQACLQQ
WIKSSDTRCCELCKYEFIMETKLKPLRKWEKLQMTSS
ERR

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 220972

Gene Symbol: MARCH8

Gene Alias: MARCH-VIII, MIR, RNF178, c-MIR