

Datasheet

ZNF485 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00220992-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF485 protein.

Immunogen: ZNF485 (NP_660355, 1 a.a. ~ 402 a.a) full-length human protein.

Sequence:

MLENYGNLVS VGLLSKPKLITQLEQGAEPWTEVREA
PSGTHAVEDYWFETKMSALKQSTSEASVLGERTKSV
MMEKGLDWEGRSSTEKNYKCEGKVFKNSSSFISH
QRNHTSEKPHKCECGIAFMNSSLLNHHKVHAGKQP
YRCIECGKFLKKHSTFINHQRIHSREKPHKCECGKTFR
KNSILLSHQRIHTGQKPYKCND CGKAFAQNAALTRHE
RIHSGEKPFKCNKCGRAFRDNSTVLEHQIHTGEKPY
QCNECGKA FRKSSTLISHQRMHTGEKPYHCSKCGKS
FRYSSSFAGHQKTHSGNKPYQCRDCGKAFTKSSTLT
GHQRIHTGEKPYHCKKCGKA FRHSSGLVEHQRLHTG
EKPYKCNECGKA FPRSSALKQHKKIHNKERAMKCS

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 220992

Gene Symbol: ZNF485

Gene Alias: -