

Datasheet

PNPLA2 purified MaxPab rabbit polyclonal antibody (D01P)

Gene Alias: 1110001C14Rik, ATGL, DESNUTRIN, DKFZp667M109, FP17548, PEDF-R, TTS-2.2, TTS2

Catalog Number: H00057104-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PNPLA2 protein.

Immunogen: PNPLA2 (NP_065109.1, 1 a.a. ~ 504 a.a) full-length human protein.

Sequence:

MFPREKTWNISFAGCGFLGVYYVGVASCLREHAPFLV
ANATHIYGASAGALTATALVTGVCLGEAGAKFIEVSKE
ARKRFLGPLHPSFNLVKIIRSFLKVLPADESHEHASGRL
GISLTRVSDGENVIISHFNSKDELIQANVCSGFIPVYCG
LIPPSLQGVRYVDGGISDNLPYELKNTITVSPFSGESD
ICPQDSSTNIHELRTNTSIQFNLRNLYRLSKALFPPEP
LVLREMCKQGYRDGLRFLQRNGLLRPNPLLALPPAR
PHGPEDKDQAVESAQAEDYSQLPGEDHILEHLPARLN
EALLEACVEPTDLLTTLNMLPVRLATAMMPYTLPLE
SALSFTIRLLEWLDPVPEDIRWMKEQTGSICQYLMRA
KRKLGRHLPSRLPEQVELRRVQSLPSVPLSCAAYREA
LPGWMRNNLSLGDALAKWEECQRQLLLGLFCTNVAF
PPEALRMRAPADPAPAPADPASPQHQLAGPAPLLSTP
APEARPVIGALGL

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 57104

Gene Symbol: PNPLA2