

Datasheet

CDX4 monoclonal antibody (M11), clone 2H7

Catalog Number: H00001046-M11

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant CDX4.

Clone Name: 2H7

Immunogen: CDX4 (NP_005184, 202 a.a. ~ 284 a.a)
full length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

RKSELAVNLGLSERQVKIWFQNRRAKERKMIKKKISQF
ENSGGSVQSDSDSISPGELPNTFFTPSAVRGFQPIEI
QQVIVSE*

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 1046

Gene Symbol: CDX4

Gene Alias: -