

## Datasheet

### ADCK5 polyclonal antibody

**Catalog Number:** PAB2648

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of ADCK5.

**Immunogen:** A synthetic peptide (conjugated with KLH) corresponding to N-terminus of human ADCK5.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** WB-Ti

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Ammonium sulfate precipitation

**Recommend Usage:** Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 203054

**Gene Symbol:** ADCK5

**Gene Alias:** FLJ35454, MGC126708

#### References:

1. Complete sequencing and characterization of 21,243 full-length human cDNAs. Ota T, Suzuki Y, Nishikawa T, Otsuki T, Sugiyama T, Irie R, Wakamatsu A, Hayashi K, Sato H, Nagai K, Kimura K, Makita H, Sekine M, Obayashi M, Nishi T, Shibahara T, Tanaka T, Ishii S,

Yamamoto J, Saito K, Kawai Y, Isono Y, Nakamura Y, Nagahari K, Murakami K, Yasuda T, Iwayanagi T, Wagatsuma M, Shiratori A, Sudo H, Hosoiri T, Kaku Y, Kodaira H, Kondo H, Sugawara M, Takahashi M, Kanda K, Yokoi T, Furuya T, Kikkawa E, Omura Y, Abe K, Kamihara K, Katsuta N, Sato K, Tanikawa M, Yamazaki M, Ninomiya K Nat Genet. 2004 Jan;36(1):40-5. Epub 2003 Dec 21.