

## Datasheet

### SRPK3 polyclonal antibody

Genomics. 1997 Aug 15;44(1):8-14.

**Catalog Number:** PAB2943**Regulatory Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of SRPK3.**Immunogen:** A synthetic peptide (conjugated with KLH) corresponding to internal region of human SRPK3.**Host:** Rabbit**Reactivity:** Human**Applications:** IHC-P, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Purification:** Protein G purification**Recommend Usage:** Western Blot (1:1000)

Immunohistochemistry (1:50-100)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 26576**Gene Symbol:** SRPK3**Gene Alias:** MGC102944, MSSK1, STK23**References:**

1. Genomic organization of two novel genes on human Xq28: compact head to head arrangement of IDH gamma and TRAP delta is conserved in rat and mouse. Brenner V, Nyakatura G, Rosenthal A, Platzer M.