

Datasheet

Progesterone polyclonal antibody

Catalog Number: PAB5159

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against progesterone.

Immunogen: 4-PREGNEN-3, 20-DIONE
3-O-CARBOXYMETHYLOXIME.

Host: Rabbit

Applications: EIA, IP

(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.