

## Datasheet

### 4-pregnen-17, 20-diol-3-one polyclonal antibody

**Catalog Number:** PAB5171

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against 4-pregnen-17, 20-diol-3-one.

**Immunogen:** 4-PREGNEN-17, 20b-DIOL-3-ONE-3-CARBOXYMETHYLOXIME.

**Host:** Rabbit

**Applications:** EIA, IP

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.