

## Datasheet

### Adamts1 polyclonal antibody

**Catalog Number:** PAB5862

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Adamts1.

**Immunogen:** A synthetic peptide corresponding to C-terminus of mouse Adamts1.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Applications:** ELISA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 11504

**Gene Symbol:** Adamts1

**Gene Alias:** ADAM-TS1, ADAMTS, ADAMTS-1, C3-C5, METH-1, METH1