

Datasheet

CD8A/CD8B monoclonal antibody, clone RFT-8

Jul;1(4):243-6.

Catalog Number: MAB5618**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody
raised against native CD8A/CD8B.**Clone Name:** RFT-8**Immunogen:** Native purified CD8A/CD8B from human
thymocytes and E-rosetted lymphocytes.**Host:** Mouse**Reactivity:** Human**Applications:** Flow Cyt, IHC-Fr
(See our web site product page for detailed applications
information)**Protocols:** See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols**Specificity:** human CD8.**Form:** Liquid**Isotype:** IgG1**Recommend Usage:** Flow Cytometry (1 ug/10⁶ cells)
The optimal working dilution should be determined by
the end user.**Storage Buffer:** In 100 mM BBS, pH 8.2**Storage Instruction:** Store at 4°C.**Entrez GeneID:** 925|926**Gene Symbol:** CD8A**Gene Alias:** CD8, Leu2, MAL, p32**References:**

1. The CD8 coreceptor revisited: one chain good, two
chains better. Zamoyska R. Immunity. 1994