

Datasheet

Dffa polyclonal antibody

Catalog Number: PAB13557

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Dffa.

Immunogen: A synthetic peptide corresponding to amino acids 314-329 of mouse Dffa.

Host: Rabbit

Reactivity: Mouse

Applications: WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes CAD.

Form: Liquid

Recommend Usage: Western Blot (1:500-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 13347

Gene Symbol: Dffa

Gene Alias: A330085O09Rik, DFF35, Dff45, ICAD, ICAD-L, ICAD-S