

Datasheet

Birc5 polyclonal antibody

Catalog Number: PAB13568

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Birc5.

Immunogen: A synthetic peptide corresponding to C-terminus of mouse Birc5.

Host: Rabbit

Reactivity: Mouse

Applications: IHC-P, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes Survivin.

Form: Liquid

Recommend Usage: Western blot (2.5-5 ug/mL)

Immunohistochemistry (10 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 11799

Gene Symbol: Birc5

Gene Alias: AAC-11, Api4, TIAP, survivin40