

Datasheet

Il17a monoclonal antibody, clone TC11-8H4 (Biotin)

Catalog Number: MAB5883

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised
against recombinant Il17a.

Clone Name: TC11-8H4

Immunogen: Recombinant protein corresponding to
mouse Il17a.

Host: Rat

Reactivity: Mouse

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: Mouse Interleukin-17 (IL-17).

Form: Liquid

Conjugation: Biotin

Isotype: IgG1, kappa

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 16171

Gene Symbol: Il17a

Gene Alias: Ctlα-8, Ctlα8, Il17