

Datasheet

PYY polyclonal antibody

Catalog Number: PAB13801

Regulation Status: For research use only (RUO)

Product Description: Guinea pig polyclonal antibody raised against synthetic peptide of PYY.

Immunogen: A synthetic peptide corresponding to pig PYY.

Host: Guinea pig

Reactivity: Cat, Mouse, Pig, Rat

Applications: IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections)

(1:1500-1:2500)

Immunohistochemistry (Frozen sections)

(1:1500-1:2500)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from antiserum

Storage Instruction: Store at 4°C on dry atmosphere.

After reconstitution with distilled water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 445018

Gene Symbol: PYY

Gene Alias: -