

Datasheet

Rotavirus NCDV polyclonal antibody

Catalog Number: PAB13817

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Rotavirus NCDV (all antigens).

Immunogen: Bovine (Nebraska Calf Diarrhea Virus).

Host: Goat

Reactivity: Viruses

Applications: IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: ICPs and late structural (virion) antigens.

Cross reactivity is >90% with human rotaviruses (reported). Uninfected cell reactivity is negative against HEp-2 cells and WI-38 cells by indirect immunofluorescence.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.