

Datasheet

Plau polyclonal antibody

Catalog Number: PAB13855

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant Plau.

Immunogen: Recombinant protein corresponding to mouse Plau.

Host: Rabbit

Reactivity: Mouse

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes high molecular weight and low molecular weight mouse Urokinase.

Form: Liquid

Recommend Usage: ELISA (1:10000)

Western Blot (1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium phosphate buffer, 100 mM NaCl, 1mM EDTA, pH 6.6

Storage Instruction: Store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 18792

Gene Symbol: Plau

Gene Alias: u-PA, uPA