

Datasheet

***Bacillus* spores polyclonal antibody**

Catalog Number: PAB13894

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against *Bacillus* spores.

Immunogen: Native purified from *Bacillus cereus* (ATCC #11778) and *Bacillus subtilis* (ATCC #9372).

Host: Rabbit

Reactivity: Bacteria

Applications: IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Reactive with spores and vegetative cells of *Bacillus cereus* and *Bacillus subtilis*. Antiserum is unabsorbed and may cross-react with other *Bacillus* species.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.