

Datasheet

Chlamydia trachomatis (Elementary Bodies) polyclonal antibody

Catalog Number: PAB13906

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against *Chlamydia trachomatis* (Elementary Bodies).

Immunogen: L2 and other serovar groups.

Host: Rabbit

Reactivity: Bacteria

Applications: ELISA, IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes purified elementary bodies, disrupted (all serovars A-K, L1-L3). Cross-reacts with *Chlamydia psittacii* and *Chlamydia pneumoniae* (TWAR). Negative uninfected cell reactivity against HEp-2 cells and egg yolk sac.

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.