

Datasheet

CML polyclonal antibody (HRP)

Catalog Number: PAB14357

Regulation Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against CML.

Immunogen: CML conjugated with KLH.

Host: Goat

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Conjugation: HRP

Recommend Usage: ELISA (1:100-1:2000)

Western Blot (1:100-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 50 mM PBS, 0.1 M NaCl, pH 7.4 and 0.01% thimerosal.

Storage Instruction: Store at -20°C on dry atmosphere.

After reconstitution with 500 uL deionized water, store at -20°C.

Aliquot to avoid repeated freezing and thawing.