

Datasheet

Apo (a) polyclonal antibody (Biotin)

Catalog Number: PAB14360

Regulation Status: For research use only (RUO)

Product Description: Sheep polyclonal antibody raised against native Apo (a).

Immunogen: Native purified human Apo (a).

Host: Sheep

Reactivity: Human

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody specifically binds to human LPA and Apo (a). No cross reaction to APOB-100/48; minor to plasminogen. This antibody can be used as second antibody in immunoblot or ELISA for detection of LPA and APOA in human plasma and lipoproteins.

Form: Lyophilized

Conjugation: Biotin

Recommend Usage: ELISA (1:100-1:10000)

Western Blot (1:100-1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 50 mM PBS, 0.10 M NaCl, 0.5 mM EDTA, 0.02% sodium azide, pH 7.4

Storage Instruction: Store at -20°C on dry atmosphere. After reconstitution with 500 uL deionized water, store at -20°C.

Aliquot to avoid repeated freezing and thawing.