

Datasheet

Cd14 polyclonal antibody

Catalog Number: PAB14382

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant Cd14.

Immunogen: Recombinant protein corresponding to mouse Cd14.

Host: Rabbit

Reactivity: Mouse

Applications: ELISA, Flow Cyt, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: The antiserum represents a marker for mouse CD14. The CD14 glycoprotein, gp 55, is present on most monocytic and macrophages like cell types: The CD14 molecule has been reported to be a receptor for endotoxin.

Form: Lyophilized

Recommend Usage: Flow Cytometry (1:3000)
The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from PBS

Storage Instruction: Store at -20°C on dry atmosphere.
After reconstitution with deionized water, store at -20°C or lower.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 12475

Gene Symbol: Cd14

Gene Alias: -