

Datasheet

DAAM2 monoclonal antibody (M01), clone 1D8

Catalog Number: H00023500-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant DAAM2.

Clone Name: 1D8

Immunogen: DAAM2 (AAH47575.1, 1 a.a. ~ 132 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAPRKRSHHGLGFLCCFGGSDIPEINLRDNHPLQFME
FSSPIPNAAEELNIRFAELVDELDLTDKNREAMFALPPEK
KWQIYCSKKKVPSLTPLATSQGSWHGVALAALACSCI
HLMFITCQPCSRCWRNNSE

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 23500

Gene Symbol: DAAM2

Gene Alias: KIAA0381, MGC90515, dJ90A20A.1