

Datasheet

DVL3 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00001857-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human DVL3 protein.

Immunogen: DVL3 (AAH32459, 1 a.a. ~ 716 a.a) full-length human protein.

Sequence:

MGETKIYHLDGQETPYLVKLPLPAERVTLADFKGVLQ
RPSYKFFFKSMDDDFGVVKEEISDDNAKLPCFNRRV
SWLVSAEGSHDPAPFCADNPSELPPPMERTGGIGD
SRPPSFHPHAGGGSQENLDNDTETDSLVSQAQRERPR
RRDGPEHATRLNGTAKGERRREPGGYDSSSTLMSSE
LETTSFDDSEDDSTSRFSSSTEQSSASRLMRRHKRR
RRKQKVSRIERSSSFSSITDSTMSLNITVTLNMEKYNF
LGISIVGQSNERGDGGIYIGSIMKGGAVAADGRIE
PGDMLLQVNEINFENMSNDDAVRVLREIVHKPGPITLTVAK
CWDPSPRGCFTLPRSEPIRPIDPAAWVSHTAAMTGTF
PAYGMSPSLSTITSTSSSITSSIPDTERLDDFHLSIHSD
MAAIVKAMASPEGLEVRDRMLLKITIPNAFIGSDVVD
WLYHNVEGFTDRREARKYASNLLKAGFIRHTVNKITFS
EQCYYIFGDLGCMANLSLHDHGGSSGASDQDTLAPL
PHPGAAPWPMAFPYQYPPPPHPYNPHPGFPELGYSY
GGGSASSQHSEGSRRSSGSNRSGSDRRKEKDPKAGD
SKSGSGSGSESDHTTRSSLRGPERRAPSERSGPAASE
HSHRSHHSLASSLRSHHHPYSGPPGVPLYGPPML
MMPPPPAAMGPPGAPPGRDLASVPPELTASRQSF
FRMAMGNPSEFFVDVM

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 1857

Gene Symbol: DVL3

Gene Alias: KIAA0208

Gene Summary: This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation. [provided by RefSeq]