

Datasheet

NOB1 polyclonal antibody

2. The ubiquitin system. Hershko A, Ciechanover A.
Annu Rev Biochem. 1998;67:425-79.

Catalog Number: PAB15596

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant NOB1.

Immunogen: Recombinant protein corresponding to *Saccharomyces cerevisiae* NOB1.

Host: Rabbit

Reactivity: Yeast

Applications: IP, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody react with *S. cerevisiae* NOB1.

Form: Liquid

Recommend Usage: Western Blot (1:400)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (50% glycerol, 1 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 854221

Gene Symbol: NOB1

Gene Alias: YOR29-07

References:

1. Nob1p is required for biogenesis of the 26S proteasome and degraded upon its maturation in *Saccharomyces cerevisiae*. Tone Y, Toh-E A. Genes Dev. 2002 Dec 15;16(24):3142-57.